

CR-NJ-168-8 Crigler Najjar Export LED lamp

EURO 5.750,00




For both patient type 1 or 2

You can use it from just born till teenager

For daily Phototherapy

8 x 21 Super LED, Spectral 445-450 nm

The intensity is easy adjustable on each panel from 0 till 100% (max. 180 MicroWatt/cm²/nm).



Height adjustable frame between 60 and 90 cm.

The height can increase with 60 cm till 150 cm with 8 aluminium separate legs



4 double independent LED panels (L = 160 cm). Total width is 60 cm.

100 - 240VAC, 50/60Hz, 400-600Watt

Delivered in wooden case.

It's easy to mount over a bed