


CR-NJ-36-4 Portable Travel LED Lamp for Crigler Najjar

EURO 3.750,00



For both patient type 1 or 2

You can use it for any age

The lamp is divided in 4 independent panels, each 63 x 28 cm
(4 x 36 Super LED)

The intensity is adjustable on each panel from 100 till 160 MicroWatt/cm²/nm.



Height adjustable frame between 60 and 90 cm.

The height can increase with 60 cm till 150 cm with 8 aluminium separate legs

Each panel has its own power supply

100 - 240VAC, 50/60Hz

NO fans (soundless)

Delivered in 2 travel bags.

It's easy to mount over a bed, without tools or screws

